



# DigiWeight



**Problem** Therapists at SRAL have no way of telling how much weight is offloaded from their patients during Gait Rehabilitation.

**Solution** A custom made electronic scale that attaches to the patients' harness and hangs overhead, measuring the offloaded weight and displaying a reading in pounds.

## Connection Points:



**Versatile** can be used on gait track and treadmill

**Easy to use** connect the two carabiners and flip the switch

**One output** Easy for the therapists' to read

**Light and portable** easily moved between systems

DTC Team 2 Section 14 Shannon Leckey, Jenny Chang, Senath Angelakudeti, Ben Vost

The contents of this poster were developed under a grant from the Department of Education, NIDRR grant number H133E130019. However, those contents do not necessarily represent the policy of the Department of Education, and you should not assume endorsement by the Federal Government.